

## Hysteroscopes

with integrated fiberoptic for illumination, autoclavable




### Compatibility

|                             |                                                                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Light port adapter</b>   | ACMI-, Wolf-, Storz / Olympus                                                                                                            |
| <b>Eyepiece</b>             | Ocular according DIN 58105 threaded                                                                                                      |
| <b>Objective</b>            | Tension-free objective composed-objective                                                                                                |
| <b>Image transmission</b>   | CAD-designed rodence system for optimal sharpness, colorsync and resolution                                                              |
| <b>Light-guiding system</b> | High quality light fibers                                                                                                                |
| <b>Optical glass</b>        | Distal and proximal Sapphire                                                                                                             |
| <b>Biocompatibility</b>     | All materials that come in contact with the body are biocompatible.<br>All outside metal parts are made of high quality stainless steel. |

| Ø      | working length | direction of view      |                        |
|--------|----------------|------------------------|------------------------|
|        |                | 0°                     | 30°                    |
| 4,0 mm | 302 mm         | <b>R300-040302-000</b> | <b>R300-040302-030</b> |

Ø 7,0 mm

Multi function valve trocar sleeve,  
metal sleeve, smooth



working length

95 mm R023-070095-001

Standard trocar



working length pyramidal

95 mm R020-070095-001




working length conical

95 mm R020-070095-002

### Diagnostic Sheaths

For Scope Ø 2,7 mm  
Outer Ø 4,0 mm



|  |                | Stopcock, fixed                                 |  |                | Stopcock, rotating     |
|--|----------------|-------------------------------------------------|--|----------------|------------------------|
|  | working length | 0° with Distance-Holder                         |  | working length |                        |
|  | <b>281 mm</b>  | <b>R210-040281-000</b>                          |  | <b>274 mm</b>  | <b>R211-040281-000</b> |
|  |                | 30° complete with Obturator and Distance-Holder |  |                |                        |
|  | <b>281 mm</b>  | <b>R210-040281-130</b>                          |  | <b>274 mm</b>  | <b>R211-040281-130</b> |
|  |                | 30° with Distance-Holder                        |  |                |                        |
|  | <b>281 mm</b>  | <b>R210-040281-030</b>                          |  | <b>274 mm</b>  | <b>R211-040281-030</b> |
|  |                | Obturator only                                  |  |                |                        |
|  |                | <b>R210-040305-00</b>                           |  |                |                        |
|  |                | Distance-Holder only                            |  |                |                        |
|  |                | <b>R210-040306-00</b>                           |  |                |                        |

### Diagnostic Sheaths

For Scope Ø 4,0 mm  
Outer Ø 5,0 mm

|  |                | Stopcock, fixed                                 |  |                | Stopcock, rotating     |
|--|----------------|-------------------------------------------------|--|----------------|------------------------|
|  | working length | 0° complete with Obturator and Distance-Holder  |  | working length |                        |
|  | <b>281 mm</b>  | <b>R210-050281-001</b>                          |  | <b>274 mm</b>  | <b>R211-050274-001</b> |
|  |                | 0° with Distance-Holder                         |  |                |                        |
|  | <b>281 mm</b>  | <b>R210-050274-000</b>                          |  | <b>274 mm</b>  | <b>R211-050274-000</b> |
|  |                | 30° complete with Obturator and Distance-Holder |  |                |                        |
|  | <b>281 mm</b>  | <b>R210-050274-130</b>                          |  | <b>274 mm</b>  | <b>R211-050274-130</b> |
|  |                | 30° with Distance-Holder                        |  |                |                        |
|  | <b>281 mm</b>  | <b>R210-050274-030</b>                          |  | <b>274 mm</b>  | <b>R211-050274-030</b> |
|  |                | Obturator only                                  |  |                |                        |
|  |                | <b>R210-050305-00</b>                           |  |                |                        |
|  |                | Distance-Holder only                            |  |                |                        |
|  |                | <b>R210-050306-00</b>                           |  |                |                        |

Operation Sheaths with 1 Working Channel

For Scope Ø 4,0 mm

Outer Ø 7,5 mm



|                |                 | Sheath with 1 rotating stopcock | Sheath with 2 rotating stopcocks |
|----------------|-----------------|---------------------------------|----------------------------------|
| working length | working channel | 30° OP - Sheath                 |                                  |
| 211 mm         | 7 CH.           | R 211-075211-030                | R 212-075211-030                 |

### Continuous-Flow-Sheaths

For Scope Ø 2,7 mm

Outer Ø 6,0

These sheaths have a separate Irrigation- and Working Channel.

The Outer Tube is detachable.



#### Sheath complete with 2 separate rotating stopcocks

| working length             | working channel | 0° Continuous-Flow-Sheath |
|----------------------------|-----------------|---------------------------|
| <b>193 mm</b>              | <b>5 CH.</b>    | <b>R214-060193-000</b>    |
| 30° Continuous-Flow-Sheath |                 |                           |
| <b>193 mm</b>              | <b>5 CH.</b>    | <b>R214-060193-030</b>    |

### Continuous-Flow-Sheaths

For Scope Ø 4,0 mm

Outer Ø 8,0 mm

#### Sheath complete with 2 separate rotating stopcocks

| working length             | working channel | 0° Continuous-Flow-Sheath |
|----------------------------|-----------------|---------------------------|
| <b>193 mm</b>              | <b>7 CH.</b>    | <b>R214-080193-000</b>    |
| 30° Continuous-Flow-Sheath |                 |                           |
| <b>193 mm</b>              | <b>7 CH.</b>    | <b>R214-080193-030</b>    |

Continuous-Flow Operating Sheaths  
for Hysteroscopes Ø 2,7 mm and Ø 2,9 mm



**Diagnostic-Continuous-Flow-Inner Sheath, 1 stopcock, rigid**

| working length | 0° - 300 mm            | Outer Ø |
|----------------|------------------------|---------|
| <b>260 mm</b>  | <b>R213-037260-000</b> | 4,0     |
|                | 30° - 300 mm           |         |
| <b>260 mm</b>  | <b>R213-037260-030</b> | 4,0     |



**Diagnostic-Continuous-Flow-Outer Sheath, 1 stopcock, rigid**

| working length | 0°, 30° - 300 mm, suitable for R213-037260-000, R213-037260-030 | Outer Ø |
|----------------|-----------------------------------------------------------------|---------|
| <b>240 mm</b>  | <b>R213-050260-000</b>                                          | 5,0     |



**Obturator**



suitable for R213-050260-000

**R213-050300-000**


### Continuous-Flow Operating Sheaths for Hysteroscopes Ø 2,7 mm and Ø 2,9 mm



| 1 Instrument Port |                 |                        |
|-------------------|-----------------|------------------------|
| working length    | working channel | 30° - 300 mm           |
| <b>228 mm</b>     | <b>5 CH.</b>    | <b>R212-038300-000</b> |



| Diagnostic-Continuous-Flow-Outer Sheath, 1 stopcock, rigid |                                            |  |
|------------------------------------------------------------|--------------------------------------------|--|
| working length                                             | 30° - 300 mm, suitable for R212-038300-000 |  |
| <b>215 mm</b>                                              | <b>R212-040300-000</b>                     |  |




Hysteroscopy - Operating Sheaths  
for Hysteroscopes Ø 2,7 mm and Ø 2,9 mm



**1 Instrument Port, 2 stopcock, rigid**

|                |                 |                        |
|----------------|-----------------|------------------------|
| working length | working channel | 30° - 300 mm           |
| <b>222 mm</b>  | <b>5 CH.</b>    | <b>R212-035248-101</b> |





### Continuous-Flow Operating Sheaths for Hysteroscopes Ø 4,0 mm



#### Diagnostic-Continuous-Flow-Inner Sheath, 1 stopcock, rigid

| working length |                        | Outer Ø |
|----------------|------------------------|---------|
| 0° - 300 mm    |                        |         |
| <b>260 mm</b>  | <b>R213-054260-000</b> | 5,3     |
| 30° - 300 mm   |                        |         |
| <b>260 mm</b>  | <b>R213-054260-030</b> | 5,3     |



#### Diagnostic-Continuous-Flow-Outer Sheath, 1 stopcock, rigid

| working length                                                   |                        | Outer Ø |
|------------------------------------------------------------------|------------------------|---------|
| 0°, 30° - 300 mm, suitable for R213-054260-000, R213-0542600-030 |                        |         |
| <b>240 mm</b>                                                    | <b>R213-065250-030</b> | 6,5     |



#### Obturator

suitable for R213-065250-030


**R213-065300-000**

Continuous-Flow Operating Sheaths  
for Hysteroscopes Ø 4,0 mm



**1 Instrument Port**

| working length | working channel | 30° - 300 mm           | Outer Ø |
|----------------|-----------------|------------------------|---------|
| <b>175 mm</b>  | <b>7 CH.</b>    | <b>R214-070000-030</b> | 5,4     |



**Diagnostic-Continuous-Flow-Outer Sheath, 1 stopcock, rigid**

| working length | 30° - 300 mm, suitable for R214-070000-030 | Outer Ø |
|----------------|--------------------------------------------|---------|
| <b>155 mm</b>  | <b>R214-070200-030</b>                     | 7,6     |



**Obturator**

suitable for R214-070200-030

**R214-070000-000**

WL 41 cm

Hysteroscopy Forceps,  
semi-rigid

Ø 5 [Char]

Ø 7 [Char]



Grasping Forceps, cup serrated

R292-017400-003



Cutting Biopsy Punch

R292-017400-002

R292-023400-002



PE-Forceps, cup, sharp

R292-017400-001

R292-023400-001



Alligator Forceps

R292-017400-011

R292-023400-011



Scissors, pointed tips

R292-017400-022

R292-023400-022



Grasping Forceps, cup serrated

R292-017400-004

R292-023400-004



Scissors, blunt tips

R292-017400-021

R292-023400-021

### Ring Applicator

with advanced handling: more safety for the user due to 3-Positions-System.

Position 0 = **Safety**

Position 1 = 1st ring

Position 2 = 2nd ring

Simplified assembly and disassembly due to **Click**-mechanism on the sliding bar.

Ø 6,0 mm

### First Incision

For use with Operating-Laparoscope



working length

Complete with Loading Device and Conus suitable for OP-Laparoscope: STANDARD

327 mm

R120-060330-001



Loading Device

R120-060000-002



Conus Piece

R120-060000-001

Ø 6,5 mm

### Second Incision

For use with a regular Trocar Sleeve



working length

Complete with Loading Device

300 mm

R120-065340-001



Loading Device

R120-060000-002



Conus Piece

R120-060000-001

**Ø 5 / 10 mm**

Myoma Drill



working length    Ø 5,0 mm

**330 mm    R121-050320-01**

Ø 10,0 mm

**330 mm    R121-100330-01**

### COHEN Intra Uterine Cannula

with spring-loaded fixation, one large and one small cone



working length    Complete

**170 mm    R122-250170-000**

Cone alone, Ø 19 mm, small

**R122-190056-000**


Cone alone, Ø 25 mm, large

**R122-250056-000**

Uterine Cannula allone

**R122-250170-001**

## Uterus Manipulator

**R122-200180-199**

Complete with 6 different cones

|                                                                                     | working length | Ø     |                                                                      |
|-------------------------------------------------------------------------------------|----------------|-------|----------------------------------------------------------------------|
|  | 45 mm          | 3 mm  | Uterus obturator, thin, short<br><b>R122-030045-000</b>              |
|  | 55 mm          | 3 mm  | Uterus obturator, thin, long<br><b>R122-030055-000</b>               |
|  | 80 mm          | 8 mm  | Uterus obturator, medium, long<br><b>R122-080080-000</b>             |
|  | 100 mm         | 10 mm | Uterus obturator, thick, extra long<br><b>R122-100100-000</b>        |
|  | 60 mm          | 25 mm | Cone, with short thorn<br><b>R122-250060-000</b>                     |
|  | 70 mm          | 25 mm | Cone, with long thorn<br><b>R122-250070-000</b>                      |
|                                                                                     |                |       | <b>Manipulator</b><br>alone, without cones<br><b>R122-200180-000</b> |

Coagulation electrodes, flexible, monopolar  
470 mm

Ball end



Ø 5 [Char]

R395-017450-001

Ø 7 [Char]

R395-023450-001

Ø 9 [Char]

R395-030450-001

Needle type



R395-017450-021

R395-023450-021

Loop



R395-023450-011

Bipolar electrodes, flexible  
360 mm



Ball end



Ø 5 [Char]

R396-000360-001

Needle type



Ø 5 [Char]

R396-000360-002



Ø 7,0 mm

### LOWSLEY retractor

for manipulation of the bladder  
for perineal prostatectomy



working length

**210 mm** R139-070210-001 straight




working length

**300 mm** R139-070300-011 curved

|                                                                        |         |
|------------------------------------------------------------------------|---------|
| OPTIKEN                                                                | 7 - 8   |
| TELESCOPES                                                             |         |
| ÓPTICAS                                                                |         |
| OTTICHE                                                                |         |
| <hr/>                                                                  |         |
| CYSTOSKOPIE                                                            | 9 - 20  |
| CYSTOSCOPY                                                             |         |
| CISTOSCOPIÁ                                                            |         |
| CISTOSCOPIA                                                            |         |
| <hr/>                                                                  |         |
| RESEKTOSKOPIE                                                          | 21 - 28 |
| RESECTOSCOPY                                                           |         |
| RESECTOSCOPIÁ                                                          |         |
| RESETTOSCOPIA                                                          |         |
| <hr/>                                                                  |         |
| URETHROTOM                                                             | 29 - 32 |
| URETHROTOME                                                            |         |
| URETRÓTOMO                                                             |         |
| URETROTOMO                                                             |         |
| <hr/>                                                                  |         |
| CYSTOSKOPIE FÜR 2,9 MM Ø OPTIK                                         | 33 - 38 |
| CYSTOSCOPY FOR 2,9 MM Ø OPTIC                                          |         |
| CISTOSCOPIÁ PARA ÓPTICA 2,9 MM Ø                                       |         |
| CISTOSCOPIA PER OTTICA 2,9 MM Ø                                        |         |
| <hr/>                                                                  |         |
| RESEKTOSKOPIE FÜR 2,9 MM Ø OPTIK                                       | 39 - 46 |
| RESECTOSCOPY FOR 2,9 MM Ø OPTIC                                        |         |
| RESECTOSCOPIÁ PARA ÓPTICA 2,9 MM Ø                                     |         |
| RESETTOSCOPIA PER OTTICA 2,9 MM Ø                                      |         |
| <hr/>                                                                  |         |
| URETHROTOMIE FÜR 2,9 MM Ø OPTIK                                        | 47 - 50 |
| URETHROTOMY FOR 2,9 MM Ø OPTIC                                         |         |
| URETRÓTOMO PARA ÓPTICA 2,9 MM Ø                                        |         |
| URETROTOMO PER OTTICA 2,9 MM Ø                                         |         |
| <hr/>                                                                  |         |
| PÄDIATRISCHE UROLOGIE / CYSTOSKOPIE UND RESEKTOSKOPIE, 9,5 & 11 CHAR Ø | 51 - 56 |
| PEDIATRIC UROLOGY / CYSTOSCOPY AND RESECTOSCOPY, 9,5 & 11 FR. Ø        |         |
| UROLOGÍA PEDIÁTRICA / CISTOSCOPIÁ Y RESECTOSCOPIÁ, 9,5 & 11 CHAR Ø     |         |
| UROLOGIA PEDIATRICA / CISTOSCOPIA E RESETTOSCOPIA, 9,5 & 11 CHAR Ø     |         |

|                                                                  |         |
|------------------------------------------------------------------|---------|
| PÄDIATRISCHE UROLOGIE / CYSTOSKOPIE UND RESEKTOSKOPIE, 13 CHAR Ø | 57 - 62 |
| PEDIATRIC UROLOGY / CYSTOSCOPY AND RESECTOSCOPY, 13 FR. Ø        |         |
| UROLOGÍA PEDIÁTRICA / CISTOSCOPIA Y RESECTOSCOPIA, 13 CHAR Ø     |         |
| UROLOGIA PEDIATRICA / CISTOSCOPIA E RESECTOSCOPIA, 13 CHAR Ø     |         |
| <hr/>                                                            |         |
| NEPHROSKOPIE                                                     | 63 - 66 |
| NEPHROSCOPY                                                      |         |
| NEFROSCOPIA                                                      |         |
| NEFROSCOPIA                                                      |         |
| <hr/>                                                            |         |
| URETERO-RENOSKOPIE                                               | 67 - 72 |
| URETERO-RENOSCOPY                                                |         |
| URETERO-RENOSCOPIA                                               |         |
| URETERO-RENOSCOPIA                                               |         |
| <hr/>                                                            |         |
| ZUBEHÖR                                                          | 73 - 78 |
| ACCESSORIES                                                      |         |
| ACCESORIOS                                                       |         |
| ACCESSORI                                                        |         |
| <hr/>                                                            |         |
| NUMMERNVERZEICHNIS                                               | 79 - 81 |
| NUMERICAL INDEX                                                  |         |
| INDICE NUMÉRICO                                                  |         |
| INDICE NUMERICO                                                  |         |
| <hr/>                                                            |         |



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**OPTIKEN**

**TELESCOPES**

**ÓPTICAS**

**OTTICHE**



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)


Cystoskope 302 mm Nutzlänge, integrierte Fiberoptik Faser für Beleuchtung, autoklavierbar  
Cystoscopes 302 mm working length, incorporated fiber optic bundle for illumination, autoclavable  
Cistoscopio 302 mm longitud, con conductor de luz de fibra, esterilizable e autoclavable  
Optica Cistoscopio lunghezza utile 302 mm a vision con conduttore di luce a fibre ottiche



Blickrichtung  
Angle of view  
Dirección visual  
Angolo di visione

Ø 4,0 mm

|     |                        |
|-----|------------------------|
| 0°  | <b>R300-040302-000</b> |
| 15° | <b>R300-040302-015</b> |
| 30° | <b>R300-040302-030</b> |
| 70° | <b>R300-040302-070</b> |



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**CYSTOSKOPIE**

**CYSTOSCOPY**


**CISTOSCOPIA**

**CISTOSCOPIA**



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Cystoskop-Schäfte, 2 Hähne fest mit LL-Adapter  
 Cystoscopy sheaths, 2 stop cocks fixed with LL-adapters  
 Camisas para cistoscopia, 2 valvas nigradas adaptador de LL  
 Camicia per cistoscopia, 2 rubinetti fissi con adattatore LL



| Ø [Char] | Obturator       | Sheath          | Sheath + Obturator | Colour Code |
|----------|-----------------|-----------------|--------------------|-------------|
| 17       | R310-170000-001 | R310-170235-001 | R311-170235-001    | Yellow      |
| 19       | R310-190000-001 | R310-190235-001 | R311-190235-001    | Green       |
| 21       | R310-210000-001 | R310-210235-001 | R311-210235-001    | Red         |
| 23       | R310-230000-001 | R310-230235-001 | R311-230235-001    | Blue        |
| 25       | R310-250000-001 | R310-250235-001 | R311-250235-001    | White       |



|                                             | Ø [Char] |                 |
|---------------------------------------------|----------|-----------------|
| Sichtobturator für Cystoskopschaft          | 17       | R310-170000-003 |
| Visual obturator for cystoscope sheath      | 19       | R310-190000-003 |
| Obturador visual para camisa de cistoscopia | 21       | R310-210000-003 |
| Otturatore di visualizzazione               | 23       | R310-230000-003 |
|                                             | 25       | R310-250000-003 |


Cystoskop-Schaft, Zentralhahn mit LL-Adapter  
Cystoscopy sheath, central valve fixed with LL-adapters  
Camisas para cistoscopio, valvula central con adaptador de LL  
Camicia cistoscopio, 2 rubinetti fissi dotato di rubinetto centrale con ingresso LL




| Ø [Char] | Obturator       | Sheath          | Sheath + Obturator | Colour Code |
|----------|-----------------|-----------------|--------------------|-------------|
| 17       | R310-170000-001 | R312-170235-001 | R313-170235-001    | Yellow      |
| 19       | R310-190000-001 | R312-190235-001 | R313-190235-001    | Green       |
| 21       | R310-210000-001 | R312-210235-001 | R313-210235-001    | Red         |
| 23       | R310-230000-001 | R312-230235-001 | R313-230235-001    | Blue        |
| 25       | R310-250000-001 | R312-250235-001 | R313-250235-001    | White       |



Cystoskopie - Optikbrücken Ø 4 mm  
Cystoscopy - Telescope bridges Ø 4 mm  
Cistoscopia - Puente para optica de Ø 4 mm  
Cistoscopio - Ponticello Ø 4 mm

**R310-040000-000**

Optikbrücke für Diagnose  
Telescope bridge for diagnostic  
Puente para optica para el diagnostico  
Ponticello diagnostico per visione diretta

**R310-040000-001**

Optikbrücke mit 1 Instrumentenkanal verschließbar  
Telescope bridge with 1 instrument channel lockable  
Puente para optica con 1 canal obturable para instrumental  
Ponticello diagnostico con 1 ingresso

**R310-040000-002**

Optikbrücke mit 2 Instrumentenkanal verschließbar  
Telescope bridge with 2 instrument channel lockable  
Puente para optica, con 2 canales obturables para instrumental  
Ponticello diagnostico con 2 ingressi

Cystoskopie - Arbeitseinsatz mit Albarran-Lenkhebel  
Cystoscopy - Working-element with Albarran deflector  
Cistoscopia - vaina de trabajo con palanca conductor tipo Albarran  
Cistoscopio - Inserto operativo con leva di ALBARRAN




1 Instrumentenkanal  
1 instrument channel  
1 canal para instrumentos  
1 canale per strumenti

W/O ratchet

**R310-040000-011**

W ratchet

**R310-040000-021**



2 Instrumentenkanäle  
2 instrument channels  
2 canales para instrumentos  
2 canali per strumenti

W/O ratchet

**R310-040000-012**

W ratchet

**R310-040000-022**



NICKELL Cystoskop Adapter für Frauen  
NICKELL cystoscope adapter for women  
NICKELL adaptador cistoscopico para mujeres  
NICKELL adattatore per donne

**R319-900-01**

Optische Instrumente für Cystoskopschäfte ab 21 Char Ø  
Optical instruments suitable for cystoscope sheath up on 21 Char Ø  
Instrumental optico para vainas cistoscópicas de 21 Char Ø  
Instrumenti ottici per camicia cistoscopica 21 Char Ø



Biopsiezange  
Biopsy forceps  
Pinza biopsia  
Pinza per biopsia

**R315-021000-01**



Biopsiestanze  
Biopsy punch  
Pinza biopsia  
Pinza tagliente per biopsia

**R315-021000-06**



Schere  
Scissor  
Tijera  
Forbici

**R315-021000-021**



Faßzange  
Grasping forceps  
Pinza de agarre  
Pinza per presa

**R315-021000-011**

Optische Instrumente für Cystoskopschäfte ab 21 Char Ø  
Optical instruments suitable for cystoscope sheath up on 21 Char Ø  
Instrumental optico para vainas cistoscópicas de 21 Char Ø  
Instrumenti ottici per camicia cistoscopica 21 Char Ø

Probeexzisionszange, optisch  
Biopsy punch, optical  
Pinza de biopsia, óptica  
Pinza per biopsia con canale per ottica



30° **R315-021000-30**


70° **R315-021000-70**

---

Adapter f. Resektoskop Schäfte  
Adaptor f. resectoscope sheath  
Adaptador para vainas resectoscópicas  
Adattatore per camicie resettoscopio



**R315-901**



**R315-902**

Flexible Instrumente 40 cm  
Flexible instruments 40 cm  
Instrumentos flexible 40 cm  
Instrumenti Flessibili 40 cm

**Nutzlänge von 34 cm auf Anfrage erhältlich**  
**Length of 34 cm available on request**  
**Largo útil de 34 cm disponible a petición**  
**Lunghezza utile 34 cm disponibile su richiesta**

|                                                                                                                                                                                                                                              | Ø 5 [Char]      | Ø 6 [Char]      | Ø 7 [Char]      | Ø 9 [Char]      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------------|
|  <p>Biopsiezange<br/>Biopsy forceps<br/>Pinza de biopsia<br/>Pinza per biopsia con<br/>doppia articolazione</p>                                             | R392-017400-001 | R392-021400-001 | R392-023400-001 | R392-030400-001 |
|  <p>Löffel gezahnt<br/>Foreign body, short<br/>Cuchara dentada<br/>Pinza per la rimozione con<br/>doppia articolazione</p>                                  | R392-017400-003 | R392-021400-003 | R392-023400-003 | R392-030400-003 |
|  <p>Biopsiezange einfachbeweglich<br/>Biopsy punch single action<br/>Pinza para biopsia articulacion<br/>simple<br/>Pinza per biopsia mono articolata</p> |                 |                 | R392-023400-004 |                 |
|  <p>Löffel oval<br/>Spoon Forceps<br/>Cuchara oval<br/>Cucchiaio ovale con<br/>denti aggressivi</p>                                                       | R392-017400-005 | R392-021400-005 | R392-023400-005 |                 |
|  <p>Fasszange<br/>Grasping forceps<br/>Pinza de agarre<br/>Pinza da presa<br/>con doppia articolazione</p>                                                |                 | R392-021400-006 |                 |                 |
|  <p>Fremdkörperzange<br/>Foreign body forceps<br/>Pinza para cuerpos extraños<br/>Pinza da presa per corpi<br/>estranei, mono articolata</p>              |                 |                 | R392-023400-009 |                 |

Flexible Instrumente 40 cm  
Flexible instruments 40 cm  
Instrumentos flexible 40 cm  
Instrumenti Flessibili 40 cm

**Nutzlänge von 34 cm auf Anfrage erhältlich**  
**Length of 34 cm available on request**  
**Largo útil de 34 cm disponible a petición**  
**Lunghezza utile 34 cm disponibile su richiesta**

|                                                                                                                                                                                                                                          | Ø 5 [Char]      | Ø 6 [Char]      | Ø 7 [Char]      | Ø 9 [Char]      |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------------|
|  <p>Fremdkörper-Zange<br/>Foreign body forceps<br/>Pinza para cuerpos extraños<br/>Pinza da presa per corpi estranei<br/>con doppia articolazione</p>   |                 |                 | R392-023400-010 |                 |
|  <p>Alligatorfaßzange<br/>Alligator grasping forceps<br/>Pinza de agarre tipo caiman<br/>Pinza da presa tipo alligator con<br/>doppia articolazione</p> | R392-017400-011 | R392-021400-011 | R392-023400-011 | R392-030400-011 |
|  <p>Steinzange<br/>Stone grasping forceps<br/>Pinza para calculus<br/>Pinza da presa per calculi con<br/>doppia articolazione</p>                     |                 |                 | R392-023400-014 |                 |
|  <p>Schere<br/>Scissor<br/>Tijera<br/>Forbici mono articolata</p>                                                                                     | R392-017400-021 | R392-021400-021 | R392-023400-021 | R392-030400-021 |
|  <p>Schere doppelt beweglich<br/>Scissor, double action<br/>Tijera, articulacion doble<br/>Forbici con doppia articolazione,<br/>lama seghettata</p>  |                 |                 | R392-023400-022 |                 |
|  <p>Hakenschere<br/>Hook scissor<br/>Tijera en forma de gancho<br/>Forbici rette uncinata mono<br/>articolate</p>                                     |                 |                 | R392-023400-023 | R392-030400-023 |

Dormia Steinfänger  
Dormia stone dislodger  
Malla para rescadar calculus  
Cestello dormia

**R394-017600-004**

4 Drähte, 5 Char Ø, 600 mm Nutzlänge  
4 wires, 5 Fr. Ø, 600 mm work.length  
Extractor de calculus de 4 alambres 5 Char Ø, 600 mm de largo util  
Cestello dormia 4 filamenti, lunghezza 600 mm

**R394-017600-006**

6 Drähte, 5 Char Ø, 600 mm Nutzlänge  
6 wires, 5 Fr. Ø, 600 mm work.length  
Extractor de calculus de 6 alambres 5 Char Ø, 600 mm de largo util  
Cestello dormia 6 filamenti, lunghezza 600 mm

**R394-000000-000**

Griff für Dormia Steinfänger  
handle for Dormia stone dislodger  
mango para Malla para rescadar calculus  
Impugnatura per Cestello dormia

Koagulationselektroden, flexibel  
Coagulation electrodes, flexible  
Electrodos para coagulación, flexible  
Elettrodo per coagulazione, flessibile

|                                                                                                                                  | Ø 3 [Char]      | Ø 5 [Char]      | Ø 7 [Char]      | Ø 9 [Char]      |
|----------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------------|
|  <p>Knopf<br/>Ball end<br/>Bola<br/>Sfera</p>   | R395-010450-001 | R395-017450-001 | R395-023450-001 | R395-030450-001 |
|  <p>Nadel<br/>Needle type<br/>Aguja<br/>Ago</p> | R395-010450-021 | R395-017450-021 | R395-023450-021 |                 |
|  <p>Schlinge<br/>Loop<br/>Asa<br/>Ansa</p>    |                 | R395-017450-011 | R395-023450-011 | R395-030450-011 |

Polypektomie-Schlinge  
Polypectomy loop  
Asa para polipectomia  
Ansa per polipectomia



Ø 5 [Char]  
R395-017400-125

Ø 7 [Char]  
R395-023400-125

Handgriff mit HF Anschluß für Diathermie Schlinge  
Handle w. HF connector for diathermy loop  
Mango con conexión de alta frecuencia para asa de diatermia  
Impugnatura per dormia con cursore ad anello dotata di attacco monopolare



R395-000000-000



Lithotripsie  
Lithotripsy  
Litotripsia  
Litotripsia

**R315-250000-39 (for 25 Char)**

Steinzange  
Stone grasping forceps  
Pinza para cálculo  
Pinza per calcoli

**R315-902**

Adapter für Steinzange  
Adapter for grasping forceps  
Adaptador para el uso de la pinza para cálculo  
Adattore con rubinetto per l'uso delle pinze per calcoli

**R315-240000-59**

Arbeitselement zur Stanz-Lithotripsie; optisch, 24 Char  
Working element for punch lithotripsy; optical, 24 Fr.  
Element de trabajo para litotripsia; óptico, 24 Char  
Strumento ottico, 24 Char di lavoro operativo per litotripsia

**R313-240235-001**


Lithotripsie Schaft mit Zentralhahn 24 Char  
Lithotripsy sheath with central valve, 24 Fr.  
Vaina para litotripsia con llave central, 24 Char  
Camicia per litotripsia con rubinetto centrale, 24 Char

**R310-240000-001**

Obturator Standard; für Lithotripsie Schaft 24 Char  
Obturator standard; for lithotripsy sheath 24 Fr.  
Obturador estándar; para vaina de litotripsiá 24 Char  
Otturatore standard; per camicia litotriptisia 24 Char

**R312-240000-001**

Optikbrücke mit 1 Instrumentenkanal 24 Char  
Telescope bridge; 1 instrument channel 24 Fr.  
Puente óptico; 1 canal para instrumentos 24 Char  
Ponte per telescopio con 1 canale per strumenti 24 Char



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**RESEKTOSKOPIE**

**RESECTOSCOPY**

**RESECTOSCOPIA**

**RESECTOSCOPIA**




[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Arbeitselement für Ø 4 mm Optik  
Working element for Ø 4 mm optic  
Elemento de trabajo para óptica Ø 4 mm  
Elemento di lavoro per l'uso con ottica Ø 4 mm

"AKTIV", schneidend durch Fingerdruck  
"ACTIVE", cutting by finger action  
"ACTIVO", cortante por la presión del dedo  
Impugnatura con movimento "ATTIVO"

**R320-040225-11**

Griff geschlossen  
Closed handle  
Mango cerrado  
Manico con anello chiuso

**R320-040225-12**

Griff offen  
Opened handle  
Mango abierto  
Manico aperto

Arbeitselement "PASSIV", schneiden durch Federzug  
Working element "PASSIVE", cutting by spring action  
Elemento de trabajo "PASSIVO", cortante por balancin  
Elemento di lavoro con movimento "PASSIVO"

**R320-040225-01**

Griff geschlossen  
Closed handle  
Mango cerrado  
Manico con anello chiuso

**R320-040225-02**

Griff offen  
Opened handle  
Mango abierto  
Manico aperto

Hochfrequenzkabel, monopolar, für Geräte von:  
High frequency cable, monopolar, for devices of:  
Cable monopolar de alta frecuencia, para equipos de:  
Cavo monopolare ad alta frequenza per apparecchi:



4 mm

Eschmann alte Gerät, Erbe T-Serie old units

**R329-040000-003** 3 mtr.

**R329-040000-005** 5 mtr.



4 mm

Aesculap, Berchtold, Martin

**R329-040000-203** 3 mtr.

**R329-040000-205** 5 mtr.



5 mm

Erbe, Storz

**R329-050000-003** 3 mtr.

**R329-050000-005** 5 mtr.



8 mm

Valleylab

**R329-080000-003** 3 mtr.


**R329-080000-005** 5 mtr.

Dauerspül Resektoskop-Schäfte  
 Resectoscope sheats, continuous flow  
 Vaina para resectoscopio con flujo continuo  
 Camicia resettore, a aspirazione continua


**Ø [Char]**

Innenschaft, ohne Spülhahn mit LL-Adapter  
 Inner Sheat, without irrigation stop cock with LL-adapter  
 Vaina interior, valvula de irrigacion sin adaptador de LL  
 Camicia interna, attacco LL senza rubinetto fisso

|    |                                                 |                        |
|----|-------------------------------------------------|------------------------|
| 24 | Drehbar<br>Rotatable<br>rotable                 | <b>R320-240000-130</b> |
|    | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-240000-030</b> |
| 27 | Drehbar<br>Rotatable<br>rotable                 | <b>R320-270000-130</b> |
|    | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-270000-030</b> |


**Ø [Char]**

Außenschaft, 2 Spülhahn fest mit LL-Adapter  
 Outer Sheat, 2 irrigation stop cock with LL-adapter  
 Vaina exterior, 2 valvula de irrigacion con adaptador de LL  
 Camicia esterna, attacco LL con 2 rubinetti fissi

|      |                                                 |                        |
|------|-------------------------------------------------|------------------------|
| 26   | Drehbar<br>Rotatable<br>rotable                 | <b>R320-260000-130</b> |
|      | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-260000-030</b> |
| 28,5 | Drehbar<br>Rotatable<br>rotable                 | <b>R320-285000-130</b> |
|      | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-285000-030</b> |

Dauerspül Resektoskop-Schäfte  
Resectoscope sheaths, continuous flow  
Vaina para resectoscopio con flujo continuo  
Camicia resettore, a aspirazione continua



Ø [Char]

Kombinierte Schäfte  
Combi sheaths  
Vainas combinadas  
Camicia con combinazione

|           |                                                 |                        |
|-----------|-------------------------------------------------|------------------------|
| 24 / 26   | Drehbar<br>Rotatable<br>rotable                 | <b>R320-260240-130</b> |
|           | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-260240-030</b> |
| 27 / 28,5 | Drehbar<br>Rotatable<br>rotable                 | <b>R320-285270-130</b> |
|           | Nicht drehbar<br>Not rotatable<br>non rotatable | <b>R320-285270-030</b> |

Resektoskopie Obturatoren  
 Resectoscopy obturators  
 Obturador para resectoscopia  
 Otturatori per resettoscopia



|                     |          |                 |
|---------------------|----------|-----------------|
| Standard Obturator  | Ø [Char] |                 |
| Standard obturators | 24       | R320-240000-001 |
| Obturador estandar  | 27       | R320-270000-001 |
| Standard otturatore |          |                 |



|                           |          |                 |
|---------------------------|----------|-----------------|
| Obturator, krümmbar       | Ø [Char] |                 |
| Obturator deflectable     | 24       | R320 240000-002 |
| Obturador flexible        | 27       | R320-270000-002 |
| Otturatore mono aticolato |          |                 |



|                                                    |          |                                            |                 |
|----------------------------------------------------|----------|--------------------------------------------|-----------------|
| Sicht-Obturatot                                    | Ø [Char] |                                            |                 |
| Visual Obturator                                   | 24       | ohne Kanal<br>without channel<br>sin canal | R320-240000-003 |
| Obturador visual                                   |          | mit Kanal<br>with channel<br>con canal     | R320-240000-005 |
| Otturatore diagnostico per visualizzazione diretta | 27       | ohne Kanal<br>without channel<br>sin canal | R320-270000-003 |
|                                                    |          | mit Kanal<br>with channel<br>con canal     | R320-270000-005 |

|                                                                                     |                                         | Ø 24 [Char]            | Ø 27 [Char]            |
|-------------------------------------------------------------------------------------|-----------------------------------------|------------------------|------------------------|
|    | Schlinge retrograde gebogen             | <b>R329-240000-001</b> | <b>R329-270000-001</b> |
|                                                                                     | Loop, retrograde curved                 | (for 0° Optic)         | (for 0° Optic)         |
|                                                                                     | Electrodo de asa, curva retrogradamente | <b>R329-240030-001</b> | <b>R329-270030-001</b> |
|                                                                                     | Ansa per taglio, angolata               | (for 30° Optic)        | (for 30° Optic)        |
|    | Schlinge gerade                         | <b>R329-240010-001</b> | <b>R329-270010-001</b> |
|                                                                                     | Loop, straight                          |                        |                        |
|                                                                                     | Asa recta                               |                        |                        |
|                                                                                     | Ansa per taglio, ad asso longitudinale  |                        |                        |
|  | Kugel                                   | <b>R329-240040-003</b> | <b>R329-270040-003</b> |
|                                                                                     | Ball                                    | (Ø 3 mm)               | (Ø 3 mm)               |
|                                                                                     | Bola                                    | <b>R329-240040-005</b> | <b>R329-270040-005</b> |
|                                                                                     | sferico                                 | (Ø 5 mm)               | (Ø 5 mm)               |
|  | Rolle                                   | <b>R329-240050-003</b> | <b>R329-270050-003</b> |
|                                                                                     | Roller barrel                           | (Ø 3 mm)               | (Ø 3 mm)               |
|                                                                                     | Rollo                                   | <b>R329-240050-005</b> | <b>R329-270050-005</b> |
|                                                                                     | Rullo                                   | (Ø 5 mm)               | (Ø 5 mm)               |
|  | Vaporisationsschlinge                   | <b>R329-240055-001</b> | <b>R329-270055-001</b> |
|                                                                                     | Vaporisation Loop                       | (for 0° Optic)         | (for 0° Optic)         |
|                                                                                     | Asa de vaporisacion                     | <b>R329-240056-001</b> |                        |
|                                                                                     | Elettrodo per vaporizzazione            | (for 30° Optic)        |                        |




Monopolare Elektroden  
 Monopolar electrodes  
 Electrodo monopolar  
 Elettrodo monopolare

|                                                                                     |                                      | <b>Ø 24 [Char]</b>     | <b>Ø 27 [Char]</b>     |
|-------------------------------------------------------------------------------------|--------------------------------------|------------------------|------------------------|
|    | Vaporisationselektrode               | <b>R329-240052-003</b> | <b>R329-270052-003</b> |
|                                                                                     | Vaporisation electrode               | (Ø 3 mm)               | (Ø 3 mm)               |
|                                                                                     | Electrodo de vaporisation            | <b>R329-240052-005</b> | <b>R329-270052-005</b> |
|                                                                                     | Elettrodo a rullo per vaporizzazione | (Ø 5 mm)               | (Ø 5 mm)               |
|    | Spike                                | <b>R329-240054-003</b> |                        |
|                                                                                     | Spike electrode                      | (Ø 3 mm)               |                        |
|                                                                                     | Electrodo de Spike                   | <b>R329-240054-005</b> |                        |
|                                                                                     | Elettrodo di Spike                   | (Ø 5 mm)               |                        |
|  | Kegel                                | <b>R329-240090-003</b> |                        |
|                                                                                     | Cone                                 |                        |                        |
|                                                                                     | Cono                                 |                        |                        |
|                                                                                     | Conico                               |                        |                        |
|  | Messer                               | <b>R329-240070-001</b> | <b>R329-270070-001</b> |
|                                                                                     | Knife                                |                        |                        |
|                                                                                     | Cuchillo                             |                        |                        |
|                                                                                     | Elettrodo di coagulazione a bisturi  |                        |                        |

### R329-990



Sterilisierhülse für Elektroden  
 Sterilization tube for electrodes  
 Embudo para esterilizacion de electrodos  
 Tubo di protezione per elettrodi



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**URETHROTOM**

**URETHROTOME**


**URETRÓTOMO**

**URETROTOMO**



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Urethrotomie  
Urethrotomy  
Uretrotomia  
Uretrotomia

**R330-210201-01**

Urethrotomschaft 21 Char Ø graduert seitlicher Kanal f. 5 Char Instrumente  
Urethrotome sheath 21 Char Ø, graduated lateral channel f. 5 Fr. Instruments  
Vaina para uretrotomo de 21 Char Ø con canal graduado lateral para instrumentos de 5 Char  
Camicia per uretrotomo 21 Char Ø

**R330-230201-001**

Zusatzschaft für Urethrotom, seitlich offen  
Supplementary sleeve, lateral open  
Vaina adicional para uretrotomo, apertura lateral  
Camicia supplementare, aperta lateralmente

**R330-230201-102**

Dauerspül Zusatzschaft für Urethrotomschaft 21 Char Ø  
Continuous Flow supplementary sleeve for urethrotom-sleeve 21 Char Ø  
Vaina adicional para irrigación continua para uretrotomo de 21 Char Ø  
Camicia supplementare a lavaggio continuo 21 Char Ø per camicia per uretrotomo



#### **R330-210000-01**

Urethrotom Standard Obturator für 21 Char Ø Schaft #R330-210201-01  
Urethrotome standard obturator for 21 Char Ø sheath #R330-210201-01  
Obturador estándar para uretrótomo de 21 Char Ø vaina #R330-210201-01  
Otturatore, per camicia per uretrotomo, 21 Char Ø #R330-210201-01




#### **R330-210000-11**

Urethrotom Standard Obturator für 21 Char Ø mit Arbeitskanal  
Urethrotome standard obturator for 21 Char Ø with instrument channel  
Obturador estándar para uretrótomo de 21 Char Ø con 1 canal para instrumentos  
Otturatore, per camicia per uretrotomo, 21 Char con 1 canale operative per strumenti

#### **R339-040000-001**

Optikbrücke  
Telescope bridge  
Puente para optica  
Ponticello per uretrotomo



Strikturmesser  
Cold knife  
Cuchillo de corte frio  
Lama uretrotomo per taglio a freddo



**R339-210000-001**  
gerade  
straight  
recto  
retto



**R339-210000-002**  
Wellenschliff  
wave shaped  
filo ondulado  
affilatura ondulata




**R339-210000-011**  
rund  
round shaped  
semicircular  
semisferico



**R339-210000-021**  
ringförmig  
ring shaped  
en forma de anillo  
ad anello



**R339-210000-031**  
sichelförmig  
sickle shaped  
en forma de hoz  
ad uncino



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**CYSTOSKOPIE**  
für 2,9 mm Ø Optik

**CYSTOSCOPY**  
for 2,9 mm Ø optic

**CISTOSCOPIA**  
para óptica 2,9 mm Ø

**CISTOSCOPIA**  
per ottica 2,9 mm Ø



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Cystoskope 302 mm Nutzlänge, integrierte Fiberoptik Faser für Beleuchtung, autoklavierbar  
Cystoscopes 302 mm working length, incorporated fiber optic bundle for illumination, autoclavable  
Cistoscopio 302 mm longitud, con conductor de luz de fibra, esterilizable e autoclavable  
Cistoscopio lunghezza utile 302 mm, con conduttore di luce a fibre ottiche sterilizzabile in autoclave



Blickrichtung  
Angle of view  
Dirección visual  
Angolo di visione

Ø 2,9 mm

|     |                        |
|-----|------------------------|
| 0°  | <b>R300-029302-000</b> |
| 15° | <b>R300-029302-015</b> |
| 30° | <b>R300-029302-030</b> |
| 70° | <b>R300-029302-070</b> |

Cystoskop-Schäfte  
Cystoscopy sheaths  
Camisas para cistoscopia  
Camicie per cistoscopia


**Ø [Char] 15**

Cystoskop-Schaft, 2 Hähne fest mit LL-Adapter  
Cystoscopy sheath, 2 stop cocks fixed with LL-adapters  
Camisa para cistoscopia, 2 valvulas con adaptador de LL  
Camicia cistoscopia, 2 rubinetti fissi incluso, 2 coni LL


|                    |                        |
|--------------------|------------------------|
| Obturator          | <b>R310-150000-001</b> |
| Sheath + Obturator | <b>R310-150235-001</b> |
| Sheath             | <b>R311-150235-001</b> |



Cystoskopie - Optikbrücken  
Cystoscopy - Telescope bridges  
Cistoscopia - Puente para optica  
Cistoscopio - Ponticello

**R310-029000-000**

Optikbrücke für Diagnose  
Telescope bridge for diagnostic  
Puente para optica para el diagnostico  
Ponticello diagnostico per visione diretta

**R310-029000-001**

Optikbrücke mit 1 Instrumentenkanal verschließbar  
Telescope bridge with 1 instrument channel lockable  
Puente para optica con 1 canal obturable para instrumental  
Ponticello diagnostico con 1 canale

**R310-029000-002**

Optikbrücke mit 2 Instrumentenkanal verschließbar  
Telescope bridge with 2 instrument channel lockable  
Puente para optica con 2 canales obturables para instrumental  
Ponticello diagnostico con 2 canali

Cystoskopie - Arbeitseinsatz mit Albarran-Lenkhebel  
Cystoscopy - Working-element with Albarran deflector  
Cistoscopia - vaina de trabajo con palanca conductor tipo Albarran  
Cistoscopio - Inserto operativo con leva di ALBARRAN



**R310-029000-011**

Arbeitseinsatz mit Albarran-Lenkhebel, 1 Instrumentenkanal  
Working-element with Albarran deflector, 1 instrument channel  
Vaina de trabajo con palanca conductora tipo Albarran, 1 canal para instrumental  
Inserto operativo con leva di ALBARRAN con 1 ingresso




**R310-029000-012**

Arbeitseinsatz mit Albarran-Lenkhebel, 2 Instrumentenkanäle  
Working-element with Albarran deflector, 2 instrument channels  
Vaina de trabajo con palanca conductora tipo Albarran, 2 canales para instrumental  
Inserto operativo con leva di ALBARRAN con 2 ingressi

Flexible Instrumente 40 cm  
Flexible instruments 40 cm  
Instrumentos flexible 40 cm  
Instrumenti Flessibili 40 cm

**Nutzlänge von 34 cm auf Anfrage erhältlich**  
**Length of 34 cm available on request**  
**Largo útil de 34 cm disponible a petición**  
**Lunghezza utile 34 cm disponibile su richiesta**

|                                                                                                                                                                                                                                             | Ø 5 [Char]      | Ø 4 [Char]      | Ø 3,5 [Char]    | Ø 3 [Char]      |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|-----------------|-----------------|
|  <p>Biopsiezeange<br/>Biopsy forceps<br/>Pinza de biopsia<br/>Pinza per biopsia con<br/>doppia articolazione</p>                                           | R392-017400-001 | R392-014400-001 | R392-012400-001 | R392-010400-001 |
|  <p>Löffel oval<br/>Spoon Forceps<br/>Cuchara oval<br/>Cucchiaio ovale con<br/>denti aggressivi</p>                                                        | R392-017400-005 |                 |                 | R392-010400-005 |
|  <p>Alligatorfaßzeange<br/>Alligator grasping forceps<br/>Pinza de agarre tipo caiman<br/>Pinza da presa tipo alligator con<br/>doppia articolazione</p> | R392-017400-011 | R392-014400-011 | R392-012400-011 | R392-010400-011 |
|  <p>Schere<br/>Scissor<br/>Tijera<br/>Forbici mono articolata</p>                                                                                        | R392-017400-021 | R392-014400-021 | R392-012400-021 | R392-010400-021 |



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**RESEKTOSKOPIE**  
für 2,9 mm Ø Optik

**RESECTOSCOPY**  
for 2,9 mm Ø optic

**RESECTOSCOPIA**  
para óptica 2,9 mm Ø

**RESECTOSCOPIA**  
per ottica 2,9 mm Ø




[www.Red-Instrumente.de](http://www.Red-Instrumente.de)

Arbeitselement für Ø 2,9 mm Optik  
Working element for Ø 2,9 mm optic  
Vaina de trabajo para óptica de Ø 2,9 mm  
Elemento di lavoro per ottico Ø 2,9 mm

"AKTIV", schneidend durch Fingerdruck  
"ACTIVE", cutting by finger action  
"ACTIVO", cortante por la presión del dedo  
impugnatura con movimento "ATTIVO"

**R320-027225-11**

Griff geschlossen  
Closed handle  
Mango cerrado  
Manico con anello chiuso

**R320-027225-12**

Griff offen  
Opened handle  
Mango abierto  
Manico aperto

Arbeitselement "PASSIV", schneiden durch Federzug  
Working element "PASSIVE", cutting by spring action  
Elemento de trabajo "PASSIVO", cortante por balancin  
Elemento di lavoro con movimento "PASSIVO"

**R320-027225-01**

Griff geschlossen  
Closed handle  
Mango cerrado  
Manico con anello chiuso

**R320-027225-02**

Griff offen  
Opened handle  
Mango abierto  
Manico aperto

Hochfrequenzkabel, monopolar, für Geräte von:  
High frequency cable, monopolar, for devices of:  
Cable monopolar de alta frecuencia, para equipos de:  
Cavo monopolare ad alta frequenza per apparecchi:



4 mm

Eschmann alte Gerät, Erbe T-Serie old units

**R329-040000-003** 3 mtr.

**R329-040000-005** 5 mtr.



4 mm

Aesculap, Berchtold, Martin

**R329-040000-203** 3 mtr.

**R329-040000-205** 5 mtr.



5 mm

Erbe, Storz

**R329-050000-003** 3 mtr.

**R329-050000-005** 5 mtr.



8 mm

Valleylab

**R329-080000-003** 3 mtr.

**R329-080000-005** 5 mtr.

Dauerspül Resektoskop-Schäfte  
 Resectoscope sheats, continuous flow  
 Vaina para resectoscopio con flujo continuo  
 Camicia resettore, a aspirazione continua



Ø [Char]

Innenschaft, 1 Spülhahn fest mit LL-Adapter  
 Inner Sheat, 1 irrigation stop cock with LL-adapter  
 Vaina interior, 1 valvula de irrigacion rigida con adaptador de LL  
 Camicia interna, con 1 Rubinetto ingresso LL


19

Drehbar  
 Rotatable  
 rotatable

**R320-190000-130**

Nicht drehbar  
 Not rotatable  
 non rotatable

**R320-190000-030**



Ø [Char]

Außenschaft, 1 Spülhahn fest mit LL-Adapter  
 Outer Sheat, 1 irrigation stop cock with LL-adapter  
 Vaina exterior, 1 valvula de irrigacion rigida con adaptador de LL  
 Camicia esterna, con 1 rubinetto ingresso LL

22

Drehbar  
 Rotatable  
 rotatable

**R320-220000-130**

Nicht drehbar  
 Not rotatable  
 non rotatable

**R320-220000-030**

Dauerspül Resektoskop-Schäfte  
Resectoscope sheaths, continuous flow  
Vaina para resectoscopio con flujo continuo  
Camicia resettore, a aspirazione continua



Ø [Char]

Kombinierte Schäfte  
Combi sheaths  
Vainas combinadas  
Camicia con combinazione

19 / 22

Drehbar  
Rotatable  
rotable

**R320-220190-130**

Nicht drehbar  
Not rotatable  
non rotatable

**R320-220190-030**



Resektoskopie Obturatoren  
Resectoscopy obturators  
Obturadores para resectoscopia  
Otturatori per resettoscopia



Standard Obturator  
Standard obturators  
Obturador estandar  
Standard otturatore

Ø [Char]  
19

R320-190000-001



Sicht-Obturatot  
Visual Obturator  
Obturador visual  
Otturatore diagnostico per visualizzazione diretta

Ø [Char]  
19


ohne Kanal  
without channel  
sin canal

R320-190000-003

Monopolare Elektroden  
 Monopolar electrodes  
 Electrodo monopolares  
 Elettrodo monopolare

|                                                                                     |                                                                                                                   |                                                                                                                 |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
|    | Schlinge retrograde gebogen<br>Loop, retrograde curved<br>Asa, curva retrogradamente<br>Ansa per taglio, angolata | <b>Ø 19 [Char]</b><br><br><b>R329-190000-001</b><br>(for 0° Optic)<br><b>R329-190030-001</b><br>(for 12° Optic) |
|    | Kugel<br>Ball<br>Bola<br>Sferico                                                                                  | <b>R329-190040-002</b><br>(Ø 2 mm)<br><br><b>R329-190040-003</b><br>(Ø 3 mm)                                    |
|  | Rolle<br>Roller barrel<br>Rollo<br>Rullo                                                                          | <b>R329-190050-003</b><br>(Ø 3 mm)                                                                              |
|  | Messer<br>Knife<br>Cuchillo<br>Elettrodo di coagulazione a bisturi                                                | <b>R329-190070-001</b>                                                                                          |





**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**URETHROTOMIE**  
für 2,9 mm Ø Optik

**URETHROTOMY**  
for 2,9 mm Ø optic

**URETRÓTOMO**  
para óptica 2,9 mm Ø

**URETROTOMO**  
per ottica 2,9 mm Ø



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Urethrotomie  
Urethrotomy  
Uretrotomia  
Uretrotomia

**R330-155201-01**

Urethrotomschaft 15,5 Char Ø graduert seitlicher Kanal f. 3 Char Instrumente  
Urethrotome sheath 15,5 Char Ø, graduated lateral channel f. 3 Fr. Instruments  
Vaina para uretrotomo de 15,5 Char Ø con canal graduado lateral para instrumentos de 3. Char  
Camicia per uretrotomo 15,5 Char Ø con scale graduate dotato di canale operativo

**R330-155201-102**

Zusatzschaft für Urethrotomschaft 15,5 Char Ø, seitlich offen  
Supplementary sleeve for urethrotom-sleeve 15,5 Char Ø, lateral open  
Vaina adicional para uretrótopo de 15,5 Char Ø, apertura lateral  
Camicia supplementare per uretrotomo 15,5 Char Ø, aperta laterale

Urethrotomie  
Urethrotomy  
Uretrotomia  
Uretrotomia




### **R330-155000-11**

Urethrotom Standard Obturator für 15,5 Char Ø Schaft R330-155201-01  
Urethrotome standard obturator for 15,5 Char Ø sheath R330-155201-01  
Obturador estándar para uretrótomo de 15,5 Char Ø vaina R330-155201-01  
Otturatore per camicia per uretrotomo 15,5 Char Ø R330-155201-01

---

### **R339-029000-001**

Optikbrücke mit 3 Char Ø Arbeitskanal  
Telescope bridge w. 3 Char Ø working channel  
Puente para optica con canal de trabajo de 3 Char Ø  
Ponticello uretrotomo per strumenti 3 Char Ø




Strikturmesser

Cold knife

Cuchillo de corte frio

Lama uretrotomo per taglio a freddo

**R339-155000-001**gerade  
straight  
recto  
retto**R339-155000-002**Wellenschliff  
wave shaped  
filo ondulado  
affilatura ondulata**R339-155000-011**rund  
round shaped  
semicircular  
semisferico**R339-155000-031**sichelförmig  
sickle shaped  
en forma de hoz  
ad uncino



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**PÄDIATRISCHE UROLOGIE**  
Cystoskopie und Resektoskopie  
9,5 & 11 Char Ø

**PEDIATRIC UROLOGY**  
Cystoscopy and Resectoscopy  
9,5 & 11 Fr. Ø

**UROLOGÍA PEDIÁTRICA**  
Cistoscopia y Resectoscopia  
9,5 & 11 Char Ø

**UROLOGIA PEDIATRICA**  
Cistoscopia e Resettoscopia  
9,5 & 11 Char Ø



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)



Optik  
Optic  
Óptica  
Ottica




Blickrichtung  
Angle of view  
Dirección visual  
Angolo di visione

Ø 2,0 mm

0°                    **R300-020188-000**

30°                   **R300-020188-030**

Cystoskop-Schäfte, Kinder  
Cystoscopy sheaths, Children  
Vainas para cistoscopia, pediátricos  
Camicie per cistoscopia, Pediatrico



**Ø [Char]**

Cystoskop-Schaft, 2 Hähne fest mit LL-Adapter  
Cystoscopy sheath, 2 stop cocks fixed with LL-adapters  
Camisa para cistoscopia, 2 valvulas con adaptador de LL  
Camicia per cistoscopia, 2 rubinetti fissi con adattatore LL

|     |                    |                        |
|-----|--------------------|------------------------|
| 9,5 | Obturator          | <b>R310-095000-001</b> |
|     | Sheath             | <b>R310-095120-001</b> |
|     | Sheath + Obturator | <b>R311-095120-001</b> |

---

|    |                    |                        |
|----|--------------------|------------------------|
| 11 | Obturator          | <b>R310-110000-001</b> |
|    | Sheath             | <b>R310-110120-001</b> |
|    | Sheath + Obturator | <b>R311-110120-001</b> |

Flexible Instrumente  
Flexible instruments  
Instrumentos flexible  
Instrumenti Flessibili

Instrumente 30 cm lang  
Instruments 30 cm long  
Instrumentos de 30 cm largo  
Instrumenti lunghezza 30 cm



Biopsiezange  
Biopsy forceps  
Pinza de biopsia  
Pinza per biopsia doppia articolazione

Ø 3 [Char]

**R392-010300-001**



Alligatorfaßzange  
Alligator grasping forceps  
Pinza de agarre  
Pinza da presa tipo alligator doppia articolazione

**R392-010300-011**



Schere  
Scissor  
Tijera  
Forbici mono articolata

**R392-010300-021**




Löffel oval  
Spoon Forceps  
Cuchara oval  
Cucchiaino ovale con  
denti aggressivi

**R392-010300-005**

Resektoskop für Kinder  
Resectoscope for children  
Resectoscopio para niños  
Resettore pediatrico

Arbeitselement  
Working element  
Vaina de trabajo  
Elemento di lavoro



**R320-027110-01**  
(passiv/passive/pasivo/passivo)



**R320-110110-030**  
Resektoskop-Schaft Ø 11 Char  
Resectoscope sheath Ø 11 Char  
Vaina para resectoscopio de Ø 11 Char  
Camicia resettore Ø 11 Char

**R320-110000-001**  
Obturator Ø 11 Char  
Obturator Ø 11 Char  
Obturador Ø 11 Char  
Otturatore Ø 11 Char

Monopolare Elektroden  
Monopolar electrodes  
Electrodos monopolares  
Elettrodo Monopolare

Ø 11 [Char]

**R329-110000-001**



Schlinge retrograde gebogen  
Loop, retrograde curved  
Asa, curva retrogradamente  
Ansa per taglio, angolata

**R329-110000-002**



Schlinge retrograde gebogen #2  
Loop, retrograde curved #2  
Asa, curva retrogradamente #2  
Ansa per taglio, angolata #2

**R329-110040-002**  
(Ø 2 mm)



Kugel  
Ball  
Bola  
Sferico

**R329-110070-001**




Messer  
Knife  
Cuchillo  
Elettrodo di coagulazione a bisturi

Urethrotomie Kleinkinder  
Urethrotomy infant  
Uretrotomia para niños  
Uretrotomia pediátrico

System für Ø 2,0 mm Optiken  
System for Ø 2,0 mm telescopes


**R330-110086-01**



11 Fr.

Uretrotom Schaft, 1 drehbarer Hahn  
Urethrotome sheath, 1 turnable stopcock  
Vaina para uretrótopo, 1 llave rotatable  
Camicia per uretrotomo, 1 rubinetto girevole

**R330-110000-01**




11 Fr.

Standard Obturator  
Obturator standard  
Obturador estándar  
Otturatore standard


**R330-020110-01**

Arbeitselement, passiv  
Working element, passive  
Vaina de trabajo, pasivo  
Elemento di lavoro, passivo



11 Fr.


**R339-110000-01**



11 Fr.

gerade / straight / recto / retto


**R339-110000-11**



11 Fr.

rund / round shaped / semicircular / semisferico


**R339-110000-41**



11 Fr.


sichelförmig / sickle shaped / en forma de hoz / ad uncino

**R339-100000-31**



11 Fr.

hakenförmig / hook shaped / en forma de gancho / ad uncino



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**PÄDIATRISCHE UROLOGIE**  
Cystoskopie und Resektoskopie  
13 Char Ø

**PEDIATRIC UROLOGY**  
Cystoscopy and Resectoscopy  
13 Fr. Ø

**UROLOGÍA PEDIÁTRICA**  
Cistoscopia y Resectoscopia  
13 Char Ø

**UROLOGIA PEDIATRICA**  
Cistoscopia e Resettoscopia  
13 Char Ø



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Optik  
Optic  
Óptica  
Ottica



Blickrichtung  
Angle of view  
Dirección visual  
Angolo di visione

Ø 2,7 mm / length 187 mm

|     |                        |
|-----|------------------------|
| 0°  | <b>R300-027188-000</b> |
| 30° | <b>R300-027188-030</b> |
| 70° | <b>R300-027188-070</b> |


Cystoskop-Schäfte, Kinder  
Cystoscopy sheaths, Children  
Vainas para cistoscopia, pediátricos  
Camicie cistoscopio, Pediatrico



Ø [Char]

Cystoskop-Schaft, 2 Hähne fest mit LL-Adapter  
Cystoscopy sheath, 2 stop cocks fixed with LL-adapters  
Camisa para cistoscopia, 2 valvulas con adaptador de LL  
Camicia per cistoscopia, 2 rubinetti fissi con adattatore LL

|    |                    |                        |
|----|--------------------|------------------------|
| 13 | Obturator          | <b>R310-130000-001</b> |
|    | Sheath             | <b>R310-130120-001</b> |
|    | Sheath + Obturator | <b>R311-130120-001</b> |



Ø [Char]

Albarran Lenkhebel  
Albarran deflector  
Palanca conductor tipo Albarran  
Leva di Albarran

|    |                        |
|----|------------------------|
| 13 | <b>R310-027000-011</b> |
|----|------------------------|



Ø [Char]

Optikbrücke  
Telescope bridge  
Puente para óptica  
Ponte diagnostico

|    |                        |
|----|------------------------|
| 13 | <b>R310-027000-000</b> |
|----|------------------------|



Flexible Instrumente  
Flexible instruments  
Instrumentos flexible  
Instrumenti Flessibili

Instrumente 30 cm lang  
Instruments 30 cm long  
Instrumentos de 30 cm largo  
Instrumenti lunghezza 30 cm



Biopsiezange  
Biopsy forceps  
Pinza de biopsia  
Pinza per biopsia doppia articolazione

Ø 3 [Char]

**R392-010300-001**



Alligatorfaßzange  
Alligator grasping forceps  
Pinza de agarre  
Pinza da presa tipo alligator doppia articolazione

**R392-010300-011**



Schere  
Scissor  
Tijera  
Forbici mono articolata


**R392-010300-021**



Löffel oval  
Spoon Forceps  
Cuchara oval  
Cucchiaio ovale con  
denti aggressivi

**R392-010300-005**

Arbeitselement  
Working element  
Vaina de trabajo  
Elemente di lavoro



**R320-027110-01**  
(passiv/passive/pasivo/passivo)



**R320-130110-030**  
Resektoskop-Schaft Ø 13 Char  
Resectoscope sheath Ø 13 Char  
Vaina para resectoscopio de Ø 13 Char  
Camicia resettore Ø 13 Char

**R320-130000-001**  
Obturator Ø 13 Char  
Obturator Ø 13 Char  
Obturador Ø 13 Char  
Otturatore Ø 13 Char

Monopolare Elektroden  
Monopolar electrodes  
Electrodos Monopolares  
Monopolar elettrodo



Schlinge retrograde gebogen  
Loop, retrograde curved  
Asa, curva retrogradamente  
Ansa per taglio, angolata

Ø 13 [Char]

**R329-130000-001**



Schlinge retrograde gebogen #2  
Loop, retrograde curved #2  
Asa, curva retrogradamente #2  
Ansa per taglio, angolata #2

**R329-130000-002**



Kugel  
Ball  
Bola  
Sferico

**R329-130040-002**  
(Ø 2 mm)




Messer  
Knife  
Cuchillo  
Elettrodo di coagulazione a bisturi

**R329-130070-001**

Urethrotomie Kinder  
Urethrotomy children  
Uretrotomia para niños  
Uretrotomia pediátrico

System für Ø 2,7 mm Optiken  
System for Ø 2,7 mm telescopes


**R330-130086-01**



13 Fr.

Uretrotom Schaft, 1 drehbarer Hahn  
Urethrotome sheath, 1 turnable stopcock  
Vaina para uretrótopo, 1 llave rotatable  
Camicia per uretrotomo, 1 rubinetto girevole


**R330-130000-01**



13 Fr.

Standard Obturator  
Obturator standard  
Obturador estándar  
Otturatore standard


**R320-027110-01**



13 Fr.

Arbeitselement, passiv  
Working element, passive  
Vaina de trabajo, pasivo  
Elemento di lavoro, passivo


**R339-130000-01**



13 Fr.

gerade / straight / recto / retto


**R339-130000-11**



13 Fr.

rund / round shaped / semicircular / semisferico


**R339-130000-41**



13 Fr.


sichelförmig / sickle shaped / en forma de hoz / ad uncino

**R339-130000-31**



13 Fr.

hakenförmig / hook shaped / en forma de gancho / ad uncino



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**NEPHROSKOPIE**

**NEPHROSCOPY**

**NEFROSCOPIA**

**NEFROSCOPIA**




[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Nephroskopie  
Nephroscopy  
Nefroscopía  
Nefroscopia

Nephroskop 7,45 mm Ø, 0° Blickrichtung, autoklavierbar  
Nephroscope 7,45 mm Ø, 0° angle of view, autoclaveable  
Nefroscopio de 7,45 mm Ø, 0° Dirección visual, autoclavable  
Nefroscopio 7,45 mm Ø, 0° Angolo di visione, autoclavabile

|        |                        |
|--------|------------------------|
| 200 mm | <b>R303-075200-000</b> |
| 250 mm | <b>R303-075250-000</b> |
| 200 mm | <b>R303-075200-090</b> |
| 250 mm | <b>R303-075250-090</b> |



Nephroscope Operationsschaft, 1 Hahn drehbar  
Nephroscopy operating sheath, 1 stop cock rotatable  
Vaina para nefroscopio, 1 valvula rotatable  
Camicia per nefroscopio, 1 rubinetto girevole

|               |                        |
|---------------|------------------------|
| 20 cm 24 Char | <b>R340-240200-002</b> |
| 25 cm 24 Char | <b>R340-240250-002</b> |
| 20 cm 26 Char | <b>R340-260200-002</b> |
| 25 cm 26 Char | <b>R340-260250-002</b> |



Obturator für Nephroskop-Schaft und Fascien Dilatator mit Kanal für Führungsdrähte  
Hollow Obturator and Fascial Dilator for Nephroscop sheath  
Obturator para vaina para nefroscopio y dilatador fascial con canal para guias  
Otturatore per camicia nefroscopio e dilatatore fasciale

|               |                        |
|---------------|------------------------|
| 20 cm 24 Char | <b>R340-240200-001</b> |
| 25 cm 24 Char | <b>R340-240250-001</b> |
| 20 cm 26 Char | <b>R340-260200-001</b> |
| 25 cm 26 Char | <b>R340-260250-001</b> |

Teleskop Bougie  
Telescope Bougie  
Dilatadores telescopios  
Dilatatore per telescopio



**Ø [Char]**

|    |                        |
|----|------------------------|
| 9  | <b>R350-090000-001</b> |
| 12 | <b>R350-120000-001</b> |
| 15 | <b>R350-150000-001</b> |
| 18 | <b>R350-180000-001</b> |
| 21 | <b>R350-210000-001</b> |
| 24 | <b>R350-240000-001</b> |
| 27 | <b>R350-270000-001</b> |
| 30 | <b>R350-300000-001</b> |

**R350-000000-199**

Set für Teleskop Bougies Ø 9 - 24 Char  
Set for telescope bougies Ø 9 - 24 Char  
Set para dilatadores telescopicos Ø 9 - 24 Char  
Set per dilatatori telescopio Ø 9 - 24 Char

**R350-000200-001** (ohne Abb.)


Führungssonde, flexibel  
Guide wire, flexible  
Guia, flexible  
Sonda guida, flessibile

**R350-000100-001** (ohne Abb.)

starr  
rigid  
rigido  
rigido

Instrumente Nephroskopie  
 Instruments nephroscopy  
 Instrumentos para Nefroscopia  
 Instrumenti per Nefroscopia

|                                                                                     |                                                                                                                              | 30 cm           | 40 cm           |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|-----------------|-----------------|
|    | Biopsiezange<br>Biopsy forceps<br>Pinza de biopsia<br>Pinza per biopsia con doppia articolazione                             | R396-030300-001 | R396-030400-001 |
|    | Löffel oval<br>Spoon Forceps<br>Cuchara oval<br>Cucchiaino ovale con denti aggressivi                                        | R396-030300-003 | R396-030400-003 |
|  | Alligatorfaßzange<br>Alligator grasping forceps<br>Pinza de agarre<br>Pinza da presa tipo alligator con doppia articolazione | R396-030300-011 | R396-030400-011 |
|  | Schere<br>Scissor<br>Tijera<br>Forbici mono articolata                                                                       | R396-030300-021 | R396-030400-021 |
|  | Hakenschere<br>Hook scissor<br>Tijera en forma de gancho<br>Forbici rette uncinata mono articolate                           | R396-030300-023 | R396-030400-023 |



**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

---

**URETERO-RENOSKOPIE**

---

**URETERO-RENOGRAPHY**

---

**URETERO-RENOSCOPIA**

---

**URETERO-RENOSCOPIA**



[www.Red-Instrumente.de](http://www.Red-Instrumente.de)



Uretero-Renoskope  
Uretero-Renoscopes  
Uretero-Renoscopio  
Uretero-Renoscopio



| Ø [Char] | Ø Arbeitskanal [Char]<br>Ø Working channel [Char]<br>Ø Canal de trabajo [Char]<br>Ø Canale di lavoro [Char] | Nutzlänge<br>Working length<br>Longitud util<br>Lunghezza operativa | Blickrichtung 8°<br>Angle of view 8°<br>Dirección visual 8°<br>Angolo di visione 8° |
|----------|-------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6 - 9    | 4                                                                                                           | 310 mm                                                              | <b>R302-060310-008</b>                                                              |
| 6 - 9    | 4                                                                                                           | 425 mm                                                              | <b>R302-060425-008</b>                                                              |
| 7,5 - 10 | 5                                                                                                           | 310 mm                                                              | <b>R302-075310-008</b>                                                              |
| 7,5 - 10 | 5                                                                                                           | 425 mm                                                              | <b>R302-075425-008</b>                                                              |



|                                                                                                                                                                                                             | Ø 3 [Char] Fr.  | Ø 3,5 [Char] Fr. | Ø 4 [Char] Fr.  | Ø 5 [Char] Fr.  |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------------|-----------------|-----------------|
|  <p>Biopsiezange<br/>Biopsy forceps<br/>Pinza para biopsia<br/>Pinza per biopsia</p>                                       | R392-010600-001 | R392-012600-001  | R392-014600-001 | R392-017600-001 |
|  <p>Löffelzange gezahnt<br/>Spoon forceps serrated<br/>Pinza de cuchara dentada<br/>Pinza per la rimozione</p>           |                 | R392-012600-008  | R392-014600-008 | R392-017600-008 |
|  <p>Löffelzange lang gezahnt<br/>Spoon forceps long serrated<br/>Pinza de cuchara, larga dentada<br/>Pinza da presa</p>  |                 |                  |                 | R392-017600-006 |
|  <p>Löffelzange kurz gezahnt<br/>Spoon forceps short serrated<br/>Pinza de cuchara, corta dentada<br/>Pinza da presa</p> |                 | R392-012600-007  |                 | R392-017600-007 |
|  <p>Faßzange gezahnt<br/>Grasping forceps serrated<br/>Pinza de agarre dentada<br/>Pinza da presa</p>                    |                 | R392-012600-010  |                 |                 |

Flexible Instrumente 60 cm

Flexible instruments 60 cm

Instrumentos flexible 60 cm

Instrumenti Flessibili 60 cm

|                                                                                     |                                                                                                                                         | Ø 3 [Char] Fr.  | Ø 3,5 [Char] Fr. | Ø 4 [Char] Fr.  | Ø 5 [Char] Fr.  |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|-----------------|------------------|-----------------|-----------------|
|    | Alligatorzange<br>Alligator forceps serrated<br>Pinza de agarre<br>Pinza da presa tipo alligator                                        | R392-010600-011 | R392-012600-011  | R392-014600-011 | R392-017600-011 |
|    | Steinstrasse lang<br>Grasping forceps Steinstrasse<br>Pinza da calculos larga<br>pinza de presa                                         |                 | R392-012600-013  |                 | R392-017600-013 |
|  | Fremdkörperzange Stein<br>Foreign Body grasper stone<br>Pinza para cuerpos extraños de calculos<br>Pinza da presa per calculi           |                 | R392-012600-014  |                 | R392-017600-014 |
|  | 3-Arm Greifer<br>3-prong grasper<br>Pinza de agarre de 3 brazos<br>Pinza tripode uncinata                                               | R392-010600-015 |                  |                 | R392-017600-015 |
|  | Rattenzahnfaßzange<br>Rat tooth grasper<br>Pinza de agarre con diente de raton<br>Pinza da presa                                        | R392-010600-016 |                  | R392-014600-016 | R392-017600-016 |
|  | Fremdkörper Greifzange gefenstert<br>Foreign body grasper fenestrated<br>Pinza de agarre cuerpo extraño<br>Pinza da presa per frammenti |                 | R392-012600-017  |                 | R392-017600-017 |
|  | Ureter-Schlitzschere<br>Uretero slitting scissors<br>tijera para uréter<br>forbici per incisione dell'uretere                           |                 | R392-012600-025  |                 | R392-017600-025 |

Dormia Steinfänger  
Dormia stone dislodger  
Malla para rescadar calculus  
Cestello dormia



**R394-010600-004**

4 Drähte 3Char/1,0 mm Ø, 600 mm Nutzlänge  
4 wires 3Fr./1,0 mm Ø, 600 mm work.length,  
Malla de 4 alambre de 3 Char/1,0 mm para extraer calculos 600 mm de largo util  
Cestello dormia 4 filamenti, lunghezza 600 mm, Ø1,0 mm

**R394-017600-004**

4 Drähte 5Char/1,7 mm Ø, 600 mm Nutzlänge  
4 wires 5Fr./1,7 mm Ø, 600 mm work.length,  
Malla de 4 alambres de 5 Char/1,7 mm para extraer calculos 600 mm de largo util  
Cestello dormia 4 filamenti, lunghezza 600 mm, Ø1,7 mm

**R394-017600-006**


6 Drähte 5Char/1,7 mm Ø, 600 mm Nutzlänge  
6 wires 5Fr./1,7 mm Ø, 600 mm work.length,  
Malla de 6 alambres de 5 Char/1,7 mm para extraer calculos 600 mm de largo util  
Cestello dormia 6 filamenti, lunghezza 600 mm, Ø1,7 mm



**R394-000000-000**

Griff für Dormia Steinfänger  
handle for Dormia stone dislodger  
mango para Malla para rescadar calculus  
Impugnatura per Cestello dormia





**Reda Instrumente GmbH**  
Gänsäcker 34 - P.O. Box 131  
D-78532 Tuttlingen

**ZUBEHÖR**

**ACCESSORIES**

**ACCESORIOS**

**ACCESSORI**



[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)

Zubehör Urologie  
Accessories Urology  
Accesorios para urologia  
Accessori Urologia

**67175-08**

STOCKMANN Penisklemme  
STOCKMANN penis seizing fcps.  
STOCKMANN pinza p. el pene  
STOCKMANN pinza serrapene

**67177-11**

STRAUSS Penisklemme  
STRAUSS penis clamp  
STRAUSS pinza p. el pene  
STRAUSS pinza serrapene

Kathetereinführungsinstrumente, Dilatationsbougies  
Catheter Introducer, Dilating Bougies  
Introductor de catéteres, Bujías dilatadoras  
Mandrini per sonde uretrali, Candelette dilatatrici



67185-35  
35 cm, 13 <sup>3</sup>/<sub>4</sub>"

GUYON  
67184-39  
39 cm, 15 <sup>1</sup>/<sub>2</sub>"



| Charr. | 28 cm (11")<br>VAN BUREN | 25 cm (10")<br>GUYON | 28 cm (11")<br>DITTEL | 30 cm (11 <sup>3</sup> / <sub>4</sub> ")<br>DITTEL |
|--------|--------------------------|----------------------|-----------------------|----------------------------------------------------|
| 8      | 67201-08                 | 67205-08             | 67206-08              | 67207-08                                           |
| 9      |                          |                      | 67206-09              | 67207-09                                           |
| 10     | 67201-10                 | 67205-10             | 67206-10              | 67207-10                                           |
| 11     |                          |                      | 67206-11              | 67207-11                                           |
| 12     | 67201-12                 | 67205-12             | 67206-12              | 67207-12                                           |
| 13     |                          |                      | 67206-13              | 67207-13                                           |
| 14     | 67201-14                 | 67205-14             | 67206-14              | 67207-14                                           |
| 15     |                          |                      | 67206-15              | 67207-15                                           |
| 16     | 67201-16                 | 67205-16             | 67206-16              | 67207-16                                           |
| 17     |                          |                      | 67206-17              | 67207-17                                           |
| 18     | 67201-18                 | 67205-18             | 67206-18              | 67207-18                                           |
| 19     |                          |                      | 67206-19              | 67207-19                                           |
| 20     | 67201-20                 | 67205-20             | 67206-20              | 67207-20                                           |
| 21     |                          |                      | 67206-21              | 67207-21                                           |
| 22     | 67201-22                 | 67205-22             | 67206-22              | 67207-22                                           |
| 23     |                          |                      | 67206-23              | 67207-23                                           |
| 24     | 67201-24                 | 67205-24             | 67206-24              | 67207-24                                           |
| 25     |                          |                      | 67206-25              | 67207-25                                           |
| 26     | 67201-26                 | 67205-26             | 67206-26              | 67207-26                                           |
| 27     |                          |                      | 67206-27              | 67207-27                                           |
| 28     | 67201-28                 | 67205-28             | 67206-28              | 67207-28                                           |
| 29     |                          |                      | 67206-29              | 67207-29                                           |
| 30     | 67201-30                 | 67205-30             | 67206-30              | 67207-30                                           |
| 32     | 67201-32                 | 67205-32             |                       |                                                    |
| 34     | 67201-34                 | 67205-34             |                       |                                                    |
| 36     | 67201-36                 | 67205-36             |                       |                                                    |
| 38     |                          | 67205-38             |                       |                                                    |
| 40     |                          | 67205-40             |                       |                                                    |



OTIS Urethrotom für Erwachsene  
OTIS urethrotome for adults  
OTIS uretrótomo para adultos  
OTIS ureterotomo per adulti



Lieferumfang  
extent of delivery

**67169-00**

OTIS Dilator  
OTIS dilatation  
OTIS instrumento de dilatación  
OTIS dilatatore



3 x Aufsätze (Kugel, Flach, Konisch)  
3 x tips (ball, flat, conical)  
juego de 3 x puntas (bola, plana, cónica)  
set di 3 x punte (a sfera, piatta e conica)

**67169-01**

2 x Messer  
2 x knives  
2 x cuchillos  
2 x bisturi

1 x Box



Blasenspritzen  
Bladder springs  
Jeringa para vejiga  
Siringa vescicale



Blasenspritze  
Bladder spring  
Jeringa para vejiga  
Siringa vescicale

Ersatzzylinder  
Spare barrels  
Cilindros de repuesto  
Cilindri di ricambio

|       |                        |                        |
|-------|------------------------|------------------------|
| 50ml  | <b>R390-000050-001</b> | <b>R390-000050-999</b> |
| 75ml  | <b>R390-000075-001</b> | <b>R390-000075-999</b> |
| 100ml | <b>R390-000100-001</b> | <b>R390-000100-999</b> |
| 150ml | <b>R390-000150-001</b> | <b>R390-000150-999</b> |

Evakuator nach ELLIK  
Evacuator acc to ELLIK  
Evacuador segun ELLIK  
Evacuatore sec. ELLIK



**R390-100000-000**  
Komplett  
complete  
completo  
completo

**R390-100000-100**  
Gummiball  
Rubber ball  
Pera de goma  
Pompetta in gomma

**R390-100000-200**  
Glasbehälter  
Glass recipient  
Recipiente de vidrio  
Recipiente in vetro

**R390-100000-300**  
Schlauch mit Verriegelungs-Adapter  
Tube with locking device  
Tubo con adattador  
Tubo con adattatore

**R390-100000-400**  
Konus mit Schlauchteil  
Cone adapter with tube  
Cono con tubo  
Cono con segmento di Tubo

Zubehör Urologie  
Accessories Urology  
Accesorios para urologia  
Accessori Urologia

Luerlock-Schlauchverbinder  
luer-lock tubing connector  
luer-lock conector para el tubo  
luer-lock connettore per tubi



**R070-000060-101**  
männlich  
male  
macho  
maschio



**R070-000060-201**  
weiblich  
female  
femenino  
femminile

Luerlock-Schlauchverbinder mit Hahn  
luer-lock tubing connector with stopcock  
luer-lock conector para el tubo con llave  
luer-lock connettore per tubi; con rubinetto



**R070-000060-001**  
männlich  
male  
macho  
maschio

**R070-000060-301**  
weiblich  
female  
femenino  
femminile

Adapter für Blasenkatheter  
Adapter for bladder catheter  
Adaptador para Catéter  
Adattore per Catetere



**R390-9001**

Dichtkappe für Instrumentenkanal  
sealing cap for instrument channel  
capuchon de cierre para canal de instrumentos  
tappino in gomma; per canale strumenti



**R210-000010-011**  
1 mm

**R210-000020-011**  
2 mm

|          |      |                 |      |                 |      |
|----------|------|-----------------|------|-----------------|------|
| 67169-00 | UR76 | 67206-21        | UR75 | R300-040302-070 | UR 8 |
| 67169-01 | UR76 | 67206-22        | UR75 | R302-060310-008 | UR68 |
| 67175-08 | UR74 | 67206-23        | UR75 | R302-060425-008 | UR68 |
| 67177-11 | UR74 | 67206-24        | UR75 | R302-075310-008 | UR68 |
| 67184-39 | UR75 | 67206-25        | UR75 | R302-075425-008 | UR68 |
| 67185-35 | UR75 | 67206-26        | UR75 | R303-075200-000 | UR64 |
| 67201-08 | UR75 | 67206-27        | UR75 | R303-075200-090 | UR64 |
| 67201-10 | UR75 | 67206-28        | UR75 | R303-075250-000 | UR64 |
| 67201-12 | UR75 | 67206-29        | UR75 | R303-075250-090 | UR64 |
| 67201-14 | UR75 | 67206-30        | UR75 | R310-027000-000 | UR59 |
| 67201-16 | UR75 | 67207-08        | UR75 | R310-027000-011 | UR59 |
| 67201-18 | UR75 | 67207-09        | UR75 | R310-029000-000 | UR36 |
| 67201-20 | UR75 | 67207-10        | UR75 | R310-029000-001 | UR36 |
| 67201-22 | UR75 | 67207-11        | UR75 | R310-029000-002 | UR36 |
| 67201-24 | UR75 | 67207-12        | UR75 | R310-029000-011 | UR37 |
| 67201-26 | UR75 | 67207-13        | UR75 | R310-029000-012 | UR37 |
| 67201-28 | UR75 | 67207-14        | UR75 | R310-040000-000 | UR12 |
| 67201-30 | UR75 | 67207-15        | UR75 | R310-040000-001 | UR12 |
| 67201-32 | UR75 | 67207-16        | UR75 | R310-040000-002 | UR12 |
| 67201-34 | UR75 | 67207-17        | UR75 | R310-040000-011 | UR13 |
| 67201-36 | UR75 | 67207-18        | UR75 | R310-040000-012 | UR13 |
| 67205-08 | UR75 | 67207-19        | UR75 | R310-040000-021 | UR13 |
| 67205-10 | UR75 | 67207-20        | UR75 | R310-040000-022 | UR13 |
| 67205-12 | UR75 | 67207-21        | UR75 | R310-095000-001 | UR53 |
| 67205-14 | UR75 | 67207-22        | UR75 | R310-095120-001 | UR53 |
| 67205-16 | UR75 | 67207-23        | UR75 | R310-110000-001 | UR53 |
| 67205-18 | UR75 | 67207-24        | UR75 | R310-110120-001 | UR53 |
| 67205-20 | UR75 | 67207-25        | UR75 | R310-130000-001 | UR59 |
| 67205-22 | UR75 | 67207-26        | UR75 | R310-130120-001 | UR59 |
| 67205-24 | UR75 | 67207-27        | UR75 | R310-150000-001 | UR35 |
| 67205-26 | UR75 | 67207-28        | UR75 | R310-150235-001 | UR35 |
| 67205-28 | UR75 | 67207-29        | UR75 | R310-170000-001 | UR10 |
| 67205-30 | UR75 | 67207-30        | UR75 | R310-170000-003 | UR10 |
| 67205-32 | UR75 | R070-000060-001 | UR78 | R310-170235-001 | UR10 |
| 67205-34 | UR75 | R070-000060-101 | UR78 | R310-190000-001 | UR10 |
| 67205-36 | UR75 | R070-000060-201 | UR78 | R310-190000-003 | UR10 |
| 67205-38 | UR75 | R070-000060-301 | UR78 | R310-190235-001 | UR10 |
| 67205-40 | UR75 | R210-000010-011 | UR78 | R310-210000-001 | UR10 |
| 67206-08 | UR75 | R210-000020-011 | UR78 | R310-210000-003 | UR10 |
| 67206-09 | UR75 | R300-020188-000 | UR52 | R310-210235-001 | UR10 |
| 67206-10 | UR75 | R300-020188-030 | UR52 | R310-230000-001 | UR10 |
| 67206-11 | UR75 | R300-027188-000 | UR58 | R310-230000-003 | UR10 |
| 67206-12 | UR75 | R300-027188-030 | UR58 | R310-230235-001 | UR10 |
| 67206-13 | UR75 | R300-027188-070 | UR58 | R310-240000-001 | UR20 |
| 67206-14 | UR75 | R300-029302-000 | UR34 | R310-250000-001 | UR10 |
| 67206-15 | UR75 | R300-029302-015 | UR34 | R310-250000-003 | UR10 |
| 67206-16 | UR75 | R300-029302-030 | UR34 | R310-250235-001 | UR10 |
| 67206-17 | UR75 | R300-029302-070 | UR34 | R311-095120-001 | UR53 |
| 67206-18 | UR75 | R300-040302-000 | UR 8 | R311-110120-001 | UR53 |
| 67206-19 | UR75 | R300-040302-015 | UR 8 | R311-130120-001 | UR59 |
| 67206-20 | UR75 | R300-040302-030 | UR 8 | R311-150235-001 | UR35 |

|                 |      |                 |      |                 |      |
|-----------------|------|-----------------|------|-----------------|------|
| R311-170235-001 | UR10 | R320-220190-130 | UR43 | R329-240054-003 | UR28 |
| R311-190235-001 | UR10 | R320-240000-001 | UR26 | R329-240054-005 | UR28 |
| R311-210235-001 | UR10 | R320-240000-002 | UR26 | R329-240055-001 | UR27 |
| R311-230235-001 | UR10 | R320-240000-003 | UR26 | R329-240056-001 | UR27 |
| R311-250235-001 | UR10 | R320-240000-005 | UR26 | R329-240070-001 | UR28 |
| R312-040000-01  | UR20 | R320-240000-030 | UR24 | R329-240090-003 | UR28 |
| R312-150000-001 | UR11 | R320-240000-130 | UR24 | R329-270000-001 | UR27 |
| R312-150235-001 | UR11 | R320-260000-030 | UR24 | R329-270010-001 | UR27 |
| R312-170235-001 | UR11 | R320-260000-130 | UR24 | R329-270030-001 | UR27 |
| R312-190235-001 | UR11 | R320-260240-030 | UR25 | R329-270040-003 | UR27 |
| R312-210235-001 | UR11 | R320-260240-130 | UR25 | R329-270040-005 | UR27 |
| R312-230235-001 | UR11 | R320-270000-001 | UR26 | R329-270050-003 | UR27 |
| R312-250235-001 | UR11 | R320-270000-002 | UR26 | R329-270050-005 | UR27 |
| R313-150235-001 | UR11 | R320-270000-003 | UR26 | R329-270055-001 | UR27 |
| R313-170235-001 | UR11 | R320-270000-005 | UR26 | R329-990        | UR28 |
| R313-190235-001 | UR11 | R320-270000-030 | UR24 | R330-110000-01  | UR56 |
| R313-210235-001 | UR11 | R320-270000-130 | UR24 | R330-110086-01  | UR56 |
| R313-230235-001 | UR11 | R320-285000-030 | UR24 | R330-130000-01  | UR62 |
| R313-240235-001 | UR20 | R320-285000-130 | UR24 | R330-130086-01  | UR62 |
| R313-250235-001 | UR11 | R320-285270-030 | UR25 | R330-155000-11  | UR49 |
| R315-021000-01  | UR14 | R320-285270-130 | UR25 | R330-155201-01  | UR48 |
| R315-021000-011 | UR14 | R329-040000-003 | UR23 | R330-155201-102 | UR48 |
| R315-021000-021 | UR14 | R329-040000-005 | UR23 | R330-210000-01  | UR31 |
| R315-021000-06  | UR14 | R329-040000-203 | UR23 | R330-210000-11  | UR31 |
| R315-021000-30  | UR15 | R329-040000-205 | UR23 | R330-210201-01  | UR30 |
| R315-021000-70  | UR15 | R329-050000-003 | UR23 | R330-230201-001 | UR30 |
| R315-240000-59  | UR20 | R329-050000-005 | UR23 | R330-230201-102 | UR30 |
| R315-250000-39  | UR20 | R329-080000-003 | UR23 | R339-029000-001 | UR49 |
| R315-901        | UR15 | R329-080000-005 | UR23 | R339-040000-001 | UR31 |
| R315-902        | UR15 | R329-110000-001 | UR55 | R339-110000-001 | UR56 |
| R319-900-01     | UR13 | R329-110000-002 | UR55 | R339-110000-011 | UR56 |
| R320-027110-01  | UR55 | R329-110040-002 | UR55 | R339-110000-031 | UR56 |
| R320-027225-01  | UR40 | R329-110070-001 | UR55 | R339-110000-041 | UR56 |
| R320-027225-02  | UR40 | R329-130000-001 | UR61 | R339-130000-001 | UR62 |
| R320-027225-11  | UR40 | R329-130000-002 | UR61 | R339-130000-011 | UR62 |
| R320-027225-12  | UR40 | R329-130040-002 | UR61 | R339-130000-031 | UR62 |
| R320-040225-01  | UR22 | R329-130070-001 | UR61 | R339-130000-041 | UR62 |
| R320-040225-02  | UR22 | R329-190000-001 | UR45 | R339-155000-001 | UR50 |
| R320-040225-11  | UR22 | R329-190030-001 | UR45 | R339-155000-002 | UR50 |
| R320-040225-12  | UR22 | R329-190040-002 | UR45 | R339-155000-011 | UR50 |
| R320-110000-001 | UR55 | R329-190050-003 | UR45 | R339-155000-031 | UR50 |
| R320-110110-030 | UR55 | R329-190070-001 | UR45 | R339-210000-001 | UR32 |
| R320-130000-001 | UR61 | R329-240000-001 | UR27 | R339-210000-002 | UR32 |
| R320-130110-030 | UR61 | R329-240010-001 | UR27 | R339-210000-011 | UR32 |
| R320-190000-001 | UR44 | R329-240030-001 | UR27 | R339-210000-021 | UR32 |
| R320-190000-003 | UR44 | R329-240040-003 | UR27 | R339-210000-031 | UR32 |
| R320-190000-030 | UR42 | R329-240040-005 | UR27 | R340-240200-001 | UR64 |
| R320-190000-130 | UR42 | R329-240050-003 | UR27 | R340-240200-002 | UR64 |
| R320-220000-030 | UR42 | R329-240050-005 | UR27 | R340-240250-001 | UR64 |
| R320-220000-130 | UR42 | R329-240052-003 | UR28 | R340-240250-002 | UR64 |
| R320-220190-030 | UR43 | R329-240052-005 | UR28 | R340-260200-001 | UR64 |

|                       |      |                       |      |                       |      |
|-----------------------|------|-----------------------|------|-----------------------|------|
| R340-260200-002 ..... | UR64 | R392-012600-013 ..... | UR70 | R394-017600-004 ..... | UR18 |
| R340-260250-001 ..... | UR64 | R392-012600-014 ..... | UR70 | R394-017600-006 ..... | UR18 |
| R340-260250-002 ..... | UR64 | R392-012600-017 ..... | UR70 | R395-000000-000 ..... | UR19 |
| R350-000000-199 ..... | UR65 | R392-012600-025 ..... | UR70 | R395-010450-001 ..... | UR19 |
| R350-000100-001 ..... | UR65 | R392-012800-025 ..... | UR70 | R395-010450-021 ..... | UR19 |
| R350-000200-001 ..... | UR65 | R392-014400-001 ..... | UR38 | R395-017400-125 ..... | UR19 |
| R350-090000-001 ..... | UR65 | R392-014400-011 ..... | UR38 | R395-017450-001 ..... | UR19 |
| R350-120000-001 ..... | UR65 | R392-014400-021 ..... | UR38 | R395-017450-011 ..... | UR19 |
| R350-150000-001 ..... | UR65 | R392-014600-001 ..... | UR69 | R395-017450-021 ..... | UR19 |
| R350-180000-001 ..... | UR65 | R392-014600-008 ..... | UR69 | R395-023400-125 ..... | UR19 |
| R350-210000-001 ..... | UR65 | R392-014600-011 ..... | UR70 | R395-023450-001 ..... | UR19 |
| R350-240000-001 ..... | UR65 | R392-014600-016 ..... | UR70 | R395-023450-011 ..... | UR19 |
| R350-270000-001 ..... | UR65 | R392-017400-001 ..... | UR16 | R395-023450-021 ..... | UR19 |
| R350-300000-001 ..... | UR65 | R392-017400-003 ..... | UR16 | R395-030450-001 ..... | UR19 |
| R390-000050-001 ..... | UR77 | R392-017400-005 ..... | UR16 | R395-030450-011 ..... | UR19 |
| R390-000050-999 ..... | UR77 | R392-017400-011 ..... | UR17 | R396-030300-001 ..... | UR66 |
| R390-000075-001 ..... | UR77 | R392-017400-021 ..... | UR17 | R396-030300-003 ..... | UR66 |
| R390-000075-999 ..... | UR77 | R392-017600-001 ..... | UR69 | R396-030300-011 ..... | UR66 |
| R390-000100-001 ..... | UR77 | R392-017600-006 ..... | UR69 | R396-030300-021 ..... | UR66 |
| R390-000100-999 ..... | UR77 | R392-017600-007 ..... | UR69 | R396-030300-023 ..... | UR66 |
| R390-000150-001 ..... | UR77 | R392-017600-008 ..... | UR69 | R396-030400-001 ..... | UR66 |
| R390-000150-999 ..... | UR77 | R392-017600-011 ..... | UR70 | R396-030400-003 ..... | UR66 |
| R390-000200-001 ..... | UR77 | R392-017600-013 ..... | UR70 | R396-030400-011 ..... | UR66 |
| R390-000200-999 ..... | UR77 | R392-017600-014 ..... | UR70 | R396-030400-021 ..... | UR66 |
| R390-100000-000 ..... | UR77 | R392-017600-015 ..... | UR70 | R396-030400-023 ..... | UR66 |
| R390-100000-001 ..... | UR77 | R392-017600-016 ..... | UR70 |                       |      |
| R390-100000-100 ..... | UR77 | R392-017600-017 ..... | UR70 |                       |      |
| R390-100000-200 ..... | UR77 | R392-017600-025 ..... | UR70 |                       |      |
| R390-100000-300 ..... | UR77 | R392-021400-001 ..... | UR16 |                       |      |
| R390-100000-400 ..... | UR77 | R392-021400-003 ..... | UR16 |                       |      |
| R390-9001 .....       | UR78 | R392-021400-005 ..... | UR16 |                       |      |
| R392-010300-001 ..... | UR54 | R392-021400-006 ..... | UR16 |                       |      |
| R392-010300-005 ..... | UR54 | R392-021400-011 ..... | UR17 |                       |      |
| R392-010300-011 ..... | UR54 | R392-021400-021 ..... | UR17 |                       |      |
| R392-010300-021 ..... | UR54 | R392-023400-001 ..... | UR16 |                       |      |
| R392-010400-001 ..... | UR38 | R392-023400-003 ..... | UR16 |                       |      |
| R392-010400-005 ..... | UR38 | R392-023400-004 ..... | UR16 |                       |      |
| R392-010400-011 ..... | UR38 | R392-023400-005 ..... | UR16 |                       |      |
| R392-010400-021 ..... | UR38 | R392-023400-009 ..... | UR16 |                       |      |
| R392-010600-001 ..... | UR69 | R392-023400-011 ..... | UR17 |                       |      |
| R392-010600-011 ..... | UR70 | R392-023400-014 ..... | UR17 |                       |      |
| R392-010600-015 ..... | UR70 | R392-023400-021 ..... | UR17 |                       |      |
| R392-010600-016 ..... | UR70 | R392-023400-022 ..... | UR17 |                       |      |
| R392-012400-001 ..... | UR38 | R392-023400-023 ..... | UR17 |                       |      |
| R392-012400-011 ..... | UR38 | R392-030400-001 ..... | UR16 |                       |      |
| R392-012400-021 ..... | UR38 | R392-030400-003 ..... | UR16 |                       |      |
| R392-012600-001 ..... | UR69 | R392-030400-011 ..... | UR17 |                       |      |
| R392-012600-007 ..... | UR69 | R392-030400-021 ..... | UR17 |                       |      |
| R392-012600-008 ..... | UR69 | R392-030400-023 ..... | UR17 |                       |      |
| R392-012600-010 ..... | UR69 | R394-000000-000 ..... | UR18 |                       |      |
| R392-012600-011 ..... | UR70 | R394-010600-004 ..... | UR71 |                       |      |

Weitere Kataloge und Prospekte  
Other cataloges and prospects

Hauptkatalog Chirurgische Instrumente  
General catalog surgical instruments

Endoskopie  
Endoscopy

Titan Instrumente  
Titanium Instruments

Container Steriltechnik  
Container Sterilization Technique

HERCULES® Drahtkneifzange  
HERCULES® Wire Cutters

Neurochirurgie  
Neuro surgery

Plastische Chirurgie  
Plastic surgery

Orthopädie  
Orthopedic

Arthroskopie  
Arthroscopy



Gänsäcker 34  
P.O. Box 131  
D-78532 Tuttlingen  
Germany



Fon +49 74 62.94 45-0  
Fax +49 74 62.94 45-20  
[info@Reda-Instrumente.de](mailto:info@Reda-Instrumente.de)  
[www.Reda-Instrumente.de](http://www.Reda-Instrumente.de)